
## RUTGERS WINLAB | Wireless Information Network Laboratory

#### Background

Batteryless, ultra-low cost RFID (radio frequency identification) tags have shown promising results for locating objects in a space without attaching devices to those objects [1]. By placing RFID tags around a room, we can estimate the location of people inside a room. However due to the interference from signals bouncing around indoors, we need a large amount of data. This summer we sought to expand this work by:

- Adding more tags. We placed 219 RFID tags across four walls along with four antennas
- Adding multiple rows of tags to increase the accuracy of localization




# Challenges

- Irregularities in the room when tagging
- Low number of tag reads for two of the four antennas (~160-180 out of 219)
- The left image illustrates the RSSI (recieve signal strength indicator) heat map when one person is in the room (largest bottom) left blob). As can be seen, there is a lot of noise compared to the original experiment on the right which we tried to replicate [1].


# Magic Room: Device-Free Localization with RFID

Sam-Fone Cheung, Aly Mustafa, Shriya Das, Xiang Meng Advisors: Aggelos Bletsas and Richard Howard



## Future Work

- Utilize neural networks to analyze the images produced for enhanced precision
- Scale up room size and number of people/test different room and tag configurations
- Remove as much noise as possible from current RSSI heat maps
- Improve and stabilize room setup
- Attain better antennas/cables for increased tag reads





#### Acknowledgements

We would like to thank Dr. Aggelos Bletsas and Dr. Rich Howard for their invaluable guidance. We would also like to thank Jennifer Shane, Ivan Seskar, and the rest of the WINLAB faculty.

#### References

[1] A. Kleniatis, A. Dimitriou and A. Bletsas, "Device-Free Localization of Multiple Humans with Passive RFID and Joint RSSI-Phase Techniques," 2024 IEEE International Conference on RFID, Cambridge, MA, USA, 2024, pp. 1-6

# WINLAB